Flo 400 :

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Connexion via wifi

Sommaire

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
The screenshot shows the Flo 400 software interface. At the top, there is a header with the logo and the text "Flo 400® AUTOMATED ULTRA FAST 3D IN SITU DETECTION". Below the header is a navigation menu with five tabs: "PROTOCOLS" (highlighted with a red border), "EXPERIMENTS" (blue border), "MAINTENANCE" (green border), "WASH CYCLE" (purple border), and "LOGOUT" (black border). Below the menu is a table with the following data:

| Protocol | User | Duration | Steps | Reagents | Date | Actions |
|---------------------|-------|----------|-------|----------|------------|--|
| Hybridation in situ | admin | 21:00:00 | 12 | 4 | 16/04/2015 | Load Edit Delete |
| Immunohistochimie | admin | 06:30:00 | 10 | 6 | 26/05/2015 | Load Edit Delete |
| WashCycle IF | admin | 00:10:00 | 9 | 5 | 10/04/2015 | Load Edit Delete |

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
The screenshot shows the FIo 400 software interface. At the top, there is a header with the logo and text: "FIo 400® AUTOMATED ULTRA FAST 3D IN SITU DETECTION". Below the header is a navigation bar with tabs: "PROTOCOLS", "EXPERIMENTS", "MAINTENANCE", "WASH CYCLE", and "LOGOUT". The "EXPERIMENTS" tab is active. Below the navigation bar is a table with the following columns: "Experiment", "Protocol", "Author", "Date", and "Actions". The table contains four rows of data:

| Experiment | Protocol | Author | Date | Actions |
|-----------------------|---------------------|-----------|------------|---|
| Hybridation In Situ | Hybridation in situ | Flogentec | 2016-01-29 | Detail Delete |
| wash | wash seeDB | Flogentec | 2016-01-15 | Detail Delete |
| seedb mouse | SeeDB 23hfluxmedium | Flogentec | 2016-01-14 | Detail Delete |
| IF WM GFPNF trit DMSO | IF WM Triton DMSO | Flogentec | 2016-01-11 | Detail Delete |

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The screenshot shows the Flogentec software interface. On the left, under 'EXPERIMENT', there are fields for Protocol (Disco 2h10), Experiment (Flogentec), Author (Flogentec), Date, and Total duration. On the right, under 'REAGENTS', there is a table:

| Reagent | Tank | Volume (mL) |
|---------|------|-------------|
| THF50% | T3 | 24 |
| THF80% | T4 | 24 |
| THF100% | T1 | 50 |
| | | 50 |
| | | 20 |

A modal dialog box titled 'Loading protocol ...' is centered on the screen. It contains a progress bar showing 10% completion, and the text 'Commands : 2/20' and 'Errors : 0'. At the bottom of the dialog, there is a 'Close' button.

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Protocol loaded

The screenshot shows the 'Protocol loaded' status. It features a blue progress bar at 100%. Below the bar, it displays 'Commands : 20/20' and 'Errors : 0'. At the bottom right, there is a 'Close' button highlighted with a red rectangular box.

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The screenshot shows the main interface with three progress bars: 'Reagents' (partially filled with purple and red), 'Temperatures' (empty), and 'Progress' (empty). Below the progress bars, there are two buttons: a green 'State' button and a red 'Stop' button.

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PROTOCOLS
EXPERIMENTS
MAINTENANCE
BUY / NEEDS
LOGOUT

STEPS PROCESSING
CONFIGURATION > PROTOCOL > EXPERIMENT

PUMP

Low : 20 %

Medium : 50 %

High : 70 %

FLOWCELLS

| | FC1 | FC2 | FC3 | FC4 |
|----|-------|-------|------|------|
| T1 | 41 °C | 40 °C | 0 °C | 0 °C |
| T2 | 0 °C | 0 °C | 0 °C | 0 °C |
| T3 | 0 °C | 0 °C | 0 °C | 0 °C |
| T4 | 0 °C | 0 °C | 0 °C | 0 °C |

Save
Cancel

REAGENTS

| Reagent | Tank | Can cycle | Color | Actions |
|--|---|--|--|---|
| PBS | T1 | Yes | | Edit Delete |
| alpha-DIG | T2 | Yes | | Edit Delete |
| DMTP | T3 | No | | Edit Delete |
| <input style="width: 80%;" type="text" value="Reagent"/> | <input style="width: 80%;" type="text" value="Tank"/> | <input style="width: 80%;" type="text"/> | <input style="width: 80%;" type="text"/> | Add reagent |

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| REAGENTS | | | | |
|----------|------|-----------|-------|--------------|
| Reagent | Tank | Can cycle | Color | Actions |
| THF50% | T3 | No | | Valid Cancel |
| THF80% | T4 | No | | Edit Delete |
| THF100% | T1 | Yes | | Edit Delete |
| THF100% | T2 | Yes | | Edit Delete |
| DCM | P1 | Yes | | Edit Delete |
| Reagent | Tank | | | Add reagent |

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| Protocol | User | Duration | Steps | Reagents | Date | Actions |
|---------------------|-------|----------|-------|----------|------------|------------------|
| Hybridation in situ | admin | 21:00:00 | 12 | 4 | 16/04/2015 | Load Edit Delete |
| Immunohistochimie | admin | 06:30:00 | 10 | 6 | 26/05/2015 | Load Edit Delete |
| WashCycle IF | admin | 00:10:00 | 9 | 5 | 10/04/2015 | Load Edit Delete |

| | | | | | |
|------|-------|--|--|------------|--------|
| Name | admin | | | 29/01/2016 | Create |
|------|-------|--|--|------------|--------|

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Duration: 02:15:08

3 → Save Cancel

| Order | Step | Duration | Reagent | Cycling | Actions |
|--------|----------|--------------|---------|---------|-------------|
| 0 | Pump Low | 00:00:00 | | | Edit Delete |
| 1 | Flow | 00:30:00 | THF50% | No | Edit Delete |
| 2 | Flow | 00:30:00 | THF80% | No | Edit Delete |
| 3 | Flow | 00:30:00 | THF100% | Yes | Edit Delete |
| 4 | Flow | 00:30:00 | THF100% | Yes | Edit Delete |
| 5 | Probe | 00:10:00 | | | Edit Delete |
| 6 | Pump OFF | 00:00:00 | | | Edit Delete |
| Action | | Hour Min Sec | Reagent | Cycling | Add step |

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